

For Research Use Only

HNRNPR Recombinant antibody, PBS Only (Capture)

Catalog Number:85078-4-PBS



Basic Information

Catalog Number:

85078-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30339

GenBank Accession Number:

BC001449

GeneID (NCBI):

10236

UNIPROT ID:

O43390

Full Name:

heterogeneous nuclear
ribonucleoprotein R

Calculated MW:

71 kDa

Purification Method:

Protein A purification

CloneNo.:

242734E1

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data