

For Research Use Only

Natriuretic peptides B/NPPB Recombinant antibody, PBS Only (Detector)

Catalog Number: 85079-1-PBS



Basic Information

Catalog Number: 85079-1-PBS	GenBank Accession Number: NM_002521.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4879	CloneNo.: 242667C7
Source: Rabbit	UNIPROT ID: P16860	
Isotype: IgG	Full Name: natriuretic peptide precursor B	
	Calculated MW: 15kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01778-2, BNP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85079-3-PBS. Detection antibody: 85079-1-PBS. Standard: Eg1029. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01778-1, BNP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85079-2-PBS. Detection antibody: 85079-1-PBS. Standard: Eg1029. Range: 0.313-40 ng/mL