

For Research Use Only

MLYCD Recombinant antibody, PBS Only (Detector)

Catalog Number:85080-2-PBS



Basic Information

Catalog Number: 85080-2-PBS	GenBank Accession Number: BC000286	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23417	CloneNo.: 242753E4
Source: Rabbit	UNIPROT ID: O95822	
Isotype: IgG	Full Name: malonyl-CoA decarboxylase	
Immunogen Catalog Number: AG7334	Calculated MW: 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01785-2, MLYCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85080-1-PBS. Detection antibody: 85080-2-PBS. Standard: Ag7334. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01785-1, MLYCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85080-3-PBS. Detection antibody: 85080-2-PBS. Standard: Ag7334. Range: 0.781-100 ng/mL.