

For Research Use Only

HTR4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85083-2-PBS



Basic Information

Catalog Number: 85083-2-PBS	GenBank Accession Number: BC074755	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3360	CloneNo.: 242771F10
Source: Rabbit	UNIPROT ID: Q13639	
Isotype: IgG	Full Name: 5-hydroxytryptamine (serotonin) receptor 4	
Immunogen Catalog Number: AG15226	Calculated MW: 388 aa, 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01788-1, HTR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85083-1-PBS. Detection antibody: 85083-2-PBS. Standard: Ag15226. Range: 0.781-100 ng/mL.