

For Research Use Only

RBPMS Recombinant antibody, PBS Only (Capture)

Catalog Number: 85086-3-PBS



Basic Information

Catalog Number: 85086-3-PBS	GenBank Accession Number: BC003608	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 11030	CloneNo.: 242730H5
Source: Rabbit	UNIPROT ID: Q93062	
Isotype: IgG	Full Name: RNA binding protein with multiple splicing	
Immunogen Catalog Number: AG7256	Calculated MW: 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01817-1, RBPMS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85086-3-PBS. Detection antibody: 85086-1-PBS. Standard: Ag7256. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01817-2, RBPMS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85086-3-PBS. Detection antibody: 85086-2-PBS. Standard: Ag7256. Range: 0.781-100 ng/mL.