

For Research Use Only

PDIA3 Recombinant antibody, PBS Only (Detector)

Catalog Number:85087-3-PBS



Basic Information

Catalog Number: 85087-3-PBS	GenBank Accession Number: BC014433	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2923	CloneNo.: 242735D9
Source: Rabbit	UNIPROT ID: P30101	
Isotype: IgG	Full Name: protein disulfide isomerase family A, member 3	
Immunogen Catalog Number: AG8741	Calculated MW: 505 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01779-2, ERp57/ERp60 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85087-1-PBS. Detection antibody: 85087-3-PBS. Standard: Ag8741. Range: 0.313-40 ng/mL.