

For Research Use Only

# Complement C2 Recombinant antibody, PBS Only (Capture)

Catalog Number:85088-1-PBS



## Basic Information

Catalog Number:

85088-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000063.6

GeneID (NCBI):

717

UNIPROT ID:

P06681-1

Full Name:

complement component 2

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

242739A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01780-2, C2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85088-1-PBS. Detection antibody: 85088-3-PBS. Standard: Eg3058. Range: 0.781-100 ng/mL.