

For Research Use Only

PRPS2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85090-1-PBS



Basic Information

Catalog Number: 85090-1-PBS	GenBank Accession Number: BC119662	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5634	CloneNo.: 242657E1
Source: Rabbit	UNIPROT ID: P11908	
Isotype: IgG	Full Name: phosphoribosyl pyrophosphate synthetase 2	
Immunogen Catalog Number: AG25785		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01794-1, PRPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85090-1-PBS. Detection antibody: 85090-2-PBS. Standard: Ag25785. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01794-2, PRPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85090-3-PBS. Detection antibody: 85090-1-PBS. Standard: Ag25785. Range: 0.781-100 ng/mL