

For Research Use Only

ACSS1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85102-3-PBS



Basic Information

Catalog Number: 85102-3-PBS	GenBank Accession Number: BC039261	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84532	CloneNo.: 242361F8
Source: Rabbit	UNIPROT ID: Q9NUB1	
Isotype: IgG	Full Name: acyl-CoA synthetase short-chain family member 1	
Immunogen Catalog Number: AG10853	Calculated MW: 689 aa, 75 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01848-2, ACSS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85102-3-PBS. Detection antibody: 85102-1-PBS. Standard: Ag10853. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP01848-1, ACSS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85102-2-PBS. Detection antibody: 85102-3-PBS. Standard: Ag10853. Range: 1.563-100 ng/mL.