

For Research Use Only

PR Recombinant antibody, PBS Only (Detector)

Catalog Number: 85112-1-PBS



Basic Information

Catalog Number: 85112-1-PBS	GenBank Accession Number: NM_000926	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5241	CloneNo.: 242686B1
Source: Rabbit	UNIPROT ID: P06401	
Isotype: IgG	Full Name: progesterone receptor	
Immunogen Catalog Number: AG22986	Calculated MW: 933 aa, 99 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01857-1, Progesterone receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85112-3-PBS. Detection antibody: 85112-1-PBS. Standard: Ag22986. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01857-2, Progesterone receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85112-2-PBS. Detection antibody: 85112-1-PBS. Standard: Ag22986. Range: 1.25-80 ng/mL.