

For Research Use Only

# CNTNAP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85116-2-PBS



## Basic Information

Catalog Number:

85116-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22844

GenBank Accession Number:

BC113373

GeneID (NCBI):

26047

UNIPROT ID:

Q9UHC6

Full Name:

contactin associated protein-like 2

Calculated MW:

1331 aa, 148 kDa

Purification Method:

Protein A purification

CloneNo.:

242391F1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01839-1, Caspr2/CNTNAP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85116-1-PBS. Detection antibody: 85116-2-PBS. Standard: Ag22844. Range: 0.781-100 ng/mL