

For Research Use Only

TPX2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85118-2-PBS



Basic Information

Catalog Number: 85118-2-PBS	GenBank Accession Number: BC020207	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22974	CloneNo.: 242528C11
Source: Rabbit	UNIPROT ID: Q9ULW0	
Isotype: IgG	Full Name: TPX2, microtubule-associated, homolog (Xenopus laevis)	
Immunogen Catalog Number: AG2334	Calculated MW: 747 aa, 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01821-1, TPX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85118-2-PBS. Detection antibody: 85118-3-PBS. Standard: Ag2334. Range: 0.625-40 ng/mL