

For Research Use Only

GPX8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85120-1-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 85120-1-PBS | GenBank Accession Number: BC029424 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 493869 | CloneNo.: 242769C6 |
| Source: Rabbit | UNIPROT ID: Q8TED1 | |
| Isotype: IgG | Full Name: glutathione peroxidase 8 (putative) | |
| Immunogen Catalog Number: AG9930 | Calculated MW: 209 aa, 24 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01832-2, GPX8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85120-2-PBS. Detection antibody: 85120-1-PBS. Standard: Ag9930. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP01832-1, GPX8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85120-3-PBS. Detection antibody: 85120-1-PBS. Standard: Ag9930. Range: 3.125-200 ng/mL.