

For Research Use Only

GYS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85122-3-PBS



Basic Information

Catalog Number: 85122-3-PBS	GenBank Accession Number: BC126310	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2998	CloneNo.: 242595C4
Source: Rabbit	UNIPROT ID: P54840	
Isotype: IgG	Full Name: glycogen synthase 2 (liver)	
Immunogen Catalog Number: AG18012	Calculated MW: 703 aa, 81 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01844-1, GYS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85122-2-PBS. Detection antibody: 85122-3-PBS. Standard: Ag18012. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP01844-2, GYS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85122-1-PBS. Detection antibody: 85122-3-PBS. Standard: Ag18012. Range: 3.125-200 ng/mL.