

For Research Use Only

RGMA Recombinant antibody, PBS Only (Capture)

Catalog Number: 85124-3-PBS



Basic Information

Catalog Number: 85124-3-PBS	GenBank Accession Number: BC015886	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 56963	CloneNo.: 242649F7
Source: Rabbit	UNIPROT ID: Q96B86	
Isotype: IgG	Full Name: RGM domain family, member A	
Immunogen Catalog Number: AG3058	Calculated MW: 450 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01837-1, RGMA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85124-3-PBS. Detection antibody: 85124-2-PBS. Standard: Ag3058. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01837-2, RGMA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85124-3-PBS. Detection antibody: 85124-1-PBS. Standard: Ag3058. Range: 0.781-100 ng/mL.