

For Research Use Only

Mouse Pon2 Recombinant antibody, PBS Only (Detector)

Catalog Number:85126-3-PBS



Basic Information

Catalog Number:

85126-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_183308.2

GeneID (NCBI):

330260

UNIPROT ID:

Q62086

Full Name:

paraoxonase 2

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242673F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01843-2, MOUSE Pon2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85126-1-PBS. Detection antibody: 85126-3-PBS. Standard: Eg2741. Range: 1.563-200 ng/mL.