

For Research Use Only

# ASCC3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85130-4-PBS



## Basic Information

<b>Catalog Number:</b> 85130-4-PBS	<b>GenBank Accession Number:</b> BC050681	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 10973	<b>CloneNo.:</b> 242755F1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8N3C0	
<b>Isotype:</b> IgG	<b>Full Name:</b> activating signal cointegrator 1 complex subunit 3	
<b>Immunogen Catalog Number:</b> AG11848	<b>Calculated MW:</b> 111 aa, 13 kDa, 251 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01828-2, ASCC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85130-4-PBS. Detection antibody: 85130-3-PBS. Standard: Ag11848. Range: 0.781-100 ng/mL.