

For Research Use Only

TIMM44 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85134-2-PBS



Basic Information

Catalog Number: 85134-2-PBS	GenBank Accession Number: BC033628	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10469	CloneNo.: 242774F7
Source: Rabbit	UNIPROT ID: O43615	
Isotype: IgG	Full Name: translocase of inner mitochondrial membrane 44 homolog (yeast)	
Immunogen Catalog Number: AG4834	Calculated MW: 452 aa, 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01826-2, TIMM44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85134-1-PBS. Detection antibody: 85134-2-PBS. Standard: Ag4834. Range: 0.313-40 ng/mL.