

For Research Use Only

EPHB4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85139-2-PBS



Basic Information

Catalog Number: 85139-2-PBS	GenBank Accession Number: BC052804	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2050	CloneNo.: 242762C4
Source: Rabbit	UNIPROT ID: P54760	
Isotype: IgG	Full Name: EPH receptor B4	
Immunogen Catalog Number: AG10033	Calculated MW: 987 aa, 108 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01825-1, EPFB4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85139-2-PBS. Detection antibody: 85139-1-PBS. Standard: Ag10033. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01825-2, EPFB4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85139-2-PBS. Detection antibody: 85139-3-PBS. Standard: Ag10033. Range: 0.781-100 ng/mL