

For Research Use Only

NAP1;NCKAP1 Recombinant antibody, PBS Only (Capture)

Catalog Number:85142-1-PBS



Basic Information

Catalog Number: 85142-1-PBS	GenBank Accession Number: BC015025	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10787	CloneNo.: 242784C1
Source: Rabbit	UNIPROT ID: Q9Y2A7	
Isotype: IgG	Full Name: NCK-associated protein 1	
Immunogen Catalog Number: AG2788	Calculated MW: 129 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01841-1, NCKAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85142-1-PBS. Detection antibody: 85142-2-PBS. Standard: Ag2788. Range: 1.563-100 ng/mL.