

For Research Use Only

# GBE1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85143-4-PBS



## Basic Information

|   |   |   |
|---|---|---|
| <b>Catalog Number:</b><br>85143-4-PBS       | <b>GenBank Accession Number:</b><br>BC012098                    | <b>Purification Method:</b><br>Protein A purification |
| <b>Concentration:</b><br>1 mg/ml            | <b>GeneID (NCBI):</b><br>2632                                   | <b>CloneNo.:</b><br>242705H1                          |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q04446                                    |   |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>glucan (1,4-alpha-), branching<br>enzyme 1 |   |
| <b>Immunogen Catalog Number:</b><br>AG14154 | <b>Calculated MW:</b><br>702 aa, 80 kDa                         |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01824-1, GBE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85143-2-PBS. Detection antibody: 85143-4-PBS. Standard: Ag14154. Range: 0.781-100 ng/mL.