

For Research Use Only

# RFC1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85148-3-PBS



## Basic Information

<b>Catalog Number:</b> 85148-3-PBS	<b>GenBank Accession Number:</b> BC051751	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1000 ug/ml	<b>GeneID (NCBI):</b> 5981	<b>CloneNo.:</b> 242818D3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P35251	
<b>Isotype:</b> IgG	<b>Full Name:</b> replication factor C (activator 1) 1, 145kDa	
<b>Immunogen Catalog Number:</b> AG13585	<b>Calculated MW:</b> 128 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01865-2, RFC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85148-4-PBS. Detection antibody: 85148-3-PBS. Standard: Ag13585. Range: 0.781-100 ng/mL.