

For Research Use Only

INPP4B Recombinant antibody, PBS Only (Detector)

Catalog Number: 85151-4-PBS



Basic Information

Catalog Number:

85151-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8676

GenBank Accession Number:

BC005273

GeneID (NCBI):

8821

UNIPROT ID:

O15327

Full Name:

inositol polyphosphate-4-phosphatase, type II, 105kDa

Calculated MW:

53 aa, 6 kDa

Purification Method:

Protein A purification

CloneNo.:

242837H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01878-1, INPP4B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85151-3-PBS. Detection antibody: 85151-4-PBS. Standard: Ag8676. Range: 0.625-80 ng/mL