

For Research Use Only

ZFP36L2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85154-1-PBS



Basic Information

Catalog Number: 85154-1-PBS	GenBank Accession Number: BC005010	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 678	CloneNo.: 242776C11
Source: Rabbit	UNIPROT ID: P47974	
Isotype: IgG	Full Name: zinc finger protein 36, C3H type-like 2	
Immunogen Catalog Number: AG13515	Calculated MW: 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01862-1, ZFP36L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85154-3-PBS. Detection antibody: 85154-1-PBS. Standard: Ag13515. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01862-2, ZFP36L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85154-2-PBS. Detection antibody: 85154-1-PBS. Standard: Ag13515. Range: 0.781-100 ng/mL