

For Research Use Only

CDS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85162-1-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 85162-1-PBS | GenBank Accession Number: BC025751 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 8760 | CloneNo.: 242787A10 |
| Source: Rabbit | UNIPROT ID: O95674 | |
| Isotype: IgG | Full Name: CDP-diacylglycerol synthase (phosphatidate cytidylyltransferase) 2 | |
| Immunogen Catalog Number: AG29228 | Calculated MW: 445 aa, 51 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01883-1, CDS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85162-3-PBS. Detection antibody: 85162-1-PBS. Standard: Ag29228. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01883-2, CDS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85162-2-PBS. Detection antibody: 85162-1-PBS. Standard: Ag29228. Range: 0.625-40 ng/mL