

For Research Use Only

# HBEGF Recombinant antibody, PBS Only (Capture)

Catalog Number: 85166-2-PBS



## Basic Information

Catalog Number:

85166-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC033097

GeneID (NCBI):

1839

UNIPROT ID:

Q99075

Full Name:

heparin-binding EGF-like growth factor

Calculated MW:

208 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

242848B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01877-2, HBEGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85166-2-PBS. Detection antibody: 85166-3-PBS. Standard: Eg3862. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01877-1, HBEGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85166-2-PBS. Detection antibody: 85166-1-PBS. Standard: Eg3862. Range: 0.781-100 ng/mL.