

For Research Use Only

# ICOSLG Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85173-1-PBS



## Basic Information

Catalog Number:

85173-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC064637

GeneID (NCBI):

23308

UNIPROT ID:

O75144

Full Name:

inducible T-cell co-stimulator ligand

Calculated MW:

33 kDa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

242241B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01890-2, ICOSLG/B7-H2/CD275 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85173-1-PBS. Detection antibody: 85173-2-PBS. Standard: Eg0008. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP01890-1, ICOSLG/B7-H2/CD275 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85173-3-PBS. Detection antibody: 85173-1-PBS. Standard: Eg0008. Range: 1.563-200 ng/mL