

For Research Use Only

TAL1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85174-2-PBS



Basic Information

Catalog Number:

85174-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003189

GeneID (NCBI):

6886

UNIPROT ID:

P17542

Full Name:

T-cell acute lymphocytic leukemia 1

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

242377B5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01902-1, Human TAL1 Recombinant Matched Antibody Pair - PBS only. 85174-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag32752. 85174-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL

The mean TAL1 concentration was determined to be 17.89 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load and 25.77 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load.