

For Research Use Only

# HUWE1 Recombinant antibody, PBS Only (Detector)

Catalog Number:85177-3-PBS



## Basic Information

|   |  |   |
|---|--|---|
| <b>Catalog Number:</b><br>85177-3-PBS       | <b>GenBank Accession Number:</b><br>BC002602               | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>10075                             | <b>CloneNo.:</b><br>242702F10                         |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q7Z6Z7                               |   |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>HECT, UBA and WWE domain containing 1 |   |
| <b>Immunogen Catalog Number:</b><br>AG13763 | <b>Calculated MW:</b><br>437 aa, 482 kDa                   |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01895-2, HUWE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85177-1-PBS. Detection antibody: 85177-3-PBS. Standard: Ag13763. Range: 0.625-80 ng/mL