

For Research Use Only

# HOXA9 Recombinant antibody, PBS Only (Detector)

Catalog Number:85178-3-PBS



## Basic Information

Catalog Number: 85178-3-PBS	GenBank Accession Number: BC010023	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3205	CloneNo.: 242733E3
Source: Rabbit	UNIPROT ID: P31269	
Isotype: IgG	Full Name: homeobox A9	
Immunogen Catalog Number: AG13350	Calculated MW: 30 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01901-2, HOXA9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85178-1-PBS. Detection antibody: 85178-3-PBS. Standard: Ag13350. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01901-1, HOXA9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85178-2-PBS. Detection antibody: 85178-3-PBS. Standard: Ag13350. Range: 0.781-100 ng/mL.