

For Research Use Only

# Mouse IL-4 Recombinant antibody, PBS Only (Capture)

Catalog Number:85179-3-PBS



## Basic Information

Catalog Number:	85179-3-PBS	GenBank Accession Number:	BC027514	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	16189	CloneNo.:	242736E4
Source:	Rabbit	UNIPROT ID:	P07750		
Isotype:	IgG	Full Name:	interleukin 4		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
mouse

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01896-2, MOUSE IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85179-3-PBS. Detection antibody: 85179-1-PBS. Standard: Eg3750. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP01896-1, MOUSE IL-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85179-3-PBS. Detection antibody: 85179-2-PBS. Standard: Eg3750. Range: 1.563-200 ng/mL