

For Research Use Only

PCBP3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85180-2-PBS



Basic Information

Catalog Number: 85180-2-PBS	GenBank Accession Number: BC012061	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54039	CloneNo.: 242791E7
Source: Rabbit	UNIPROT ID: P57721	
Isotype: IgG	Full Name: poly(rC) binding protein 3	
Immunogen Catalog Number: AG9202	Calculated MW: 313 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01894-1, PCBP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85180-2-PBS. Detection antibody: 85180-3-PBS. Standard: Ag9202. Range: 0.781-100 ng/mL.