

For Research Use Only

# SHC4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85181-2-PBS



## Basic Information

<b>Catalog Number:</b> 85181-2-PBS	<b>GenBank Accession Number:</b> BC033907	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 399694	<b>CloneNo.:</b> 242805E1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q6S5L8	
<b>Isotype:</b> IgG	<b>Full Name:</b> SHC (Src homology 2 domain containing) family, member 4	
<b>Immunogen Catalog Number:</b> AG3336	<b>Calculated MW:</b> 630 aa, 69 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01891-1, SHC4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85181-2-PBS. Detection antibody: 85181-3-PBS. Standard: Ag3336. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01891-2, SHC4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85181-2-PBS. Detection antibody: 85181-1-PBS. Standard: Ag3336. Range: 0.781-100 ng/mL.