

For Research Use Only

SHC4 Recombinant antibody

Catalog Number:85181-4-RR



Basic Information

Catalog Number:

85181-4-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3336

GenBank Accession Number:

BC033907

GeneID (NCBI):

399694

UNIPROT ID:

Q6S5L8

Full Name:

SHC (Src homology 2 domain containing) family, member 4

Calculated MW:

630 aa, 69 kDa

Observed MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

242805E3

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : fetal human brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 85181-4-RR (SHC4 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours.