

For Research Use Only

MAFF Recombinant antibody, PBS Only (Capture)

Catalog Number: 85186-1-PBS



Basic Information

Catalog Number: 85186-1-PBS	GenBank Accession Number: BC015037	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23764	CloneNo.: 242859B8
Source: Rabbit	UNIPROT ID: Q9ULX9	
Isotype: IgG	Full Name: v-maf musculoaponeurotic fibrosarcoma oncogene homolog F (avian)	
Immunogen Catalog Number: AG3460	Calculated MW: 17 aa, 8 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01898-1, MAFF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85186-1-PBS. Detection antibody: 85186-2-PBS. Standard: Ag3460. Range: 0.781-100 ng/mL.