

For Research Use Only

NDUFA5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85187-2-PBS



Basic Information

Catalog Number: 85187-2-PBS	GenBank Accession Number: BC000813	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4698	CloneNo.: 242860B3
Source: Rabbit	UNIPROT ID: Q16718	
Isotype: IgG	Full Name: NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5, 13kDa	
Immunogen Catalog Number: AG9995	Calculated MW: 13 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01899-2, NDUFA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85187-3-PBS. Detection antibody: 85187-2-PBS. Standard: Ag9995. Range: 0.781-100 ng/mL.