

For Research Use Only

Caveolin-3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85190-2-PBS



Basic Information

Catalog Number: 85190-2-PBS	GenBank Accession Number: BC069368	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 859	CloneNo.: 242897B7
Source: Rabbit	UNIPROT ID: P56539	
Isotype: IgG	Full Name: caveolin 3	
Immunogen Catalog Number: AG28009	Calculated MW: 151 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01900-2, Caveolin-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85190-2-PBS. Detection antibody: 85190-1-PBS. Standard: Ag28009. Range: 0.781-100 ng/mL.