

For Research Use Only

SLC39A8 Recombinant antibody, PBS Only (Detector)

Catalog Number:85215-2-PBS



Basic Information

Catalog Number:

85215-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14815

GenBank Accession Number:

BC012125

GeneID (NCBI):

64116

UNIPROT ID:

Q9C0K1

Full Name:

solute carrier family 39 (zinc transporter), member 8

Calculated MW:

460 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

242773D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01908-1, ZIP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85215-1-PBS. Detection antibody: 85215-2-PBS. Standard: Ag14815. Range: 1.563-100 ng/mL.