

For Research Use Only

# VAMP8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85400-3-PBS



## Basic Information

Catalog Number:

85400-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7903

GenBank Accession Number:

BC001634

GeneID (NCBI):

8673

UNIPROT ID:

Q9BV40

Full Name:

vesicle-associated membrane protein

8 (endobrevin)

Calculated MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

242830G9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01931-1, VAMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85400-1-PBS. Detection antibody: 85400-3-PBS. Standard: Ag7903. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 85400-3-PBS against Human VAMP8 was performed. The affinity constant is 0.753 nM.