

For Research Use Only

DNA ligase III/LIG3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85409-2-PBS



Basic Information

Catalog Number:	85409-2-PBS	GenBank Accession Number:	BC068005	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	3980	CloneNo.:	242870E2
Source:	Rabbit	UNIPROT ID:	P49916		
Isotype:	IgG	Full Name:	ligase III, DNA, ATP-dependent		
Immunogen Catalog Number:	AG24917	Calculated MW:	100-110 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01953-1, DNA ligase III/LIG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85409-2-PBS. Detection antibody: 85409-1-PBS. Standard: Ag24917. Range: 0.781-100 ng/mL