

For Research Use Only

CHGA Recombinant antibody, PBS Only (Capture)

Catalog Number: 85605-2-PBS



Basic Information

Catalog Number: 85605-2-PBS	GenBank Accession Number: BC006459	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 1113	CloneNo.: 243048H7
Source: Rabbit	UNIPROT ID: P10645	
Isotype: IgG	Full Name: chromogranin A (parathyroid secretory protein 1)	
Immunogen Catalog Number: AG19780	Calculated MW: 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01975-2, Chromogranin A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85605-2-PBS. Detection antibody: 85605-1-PBS. Standard: Ag19780. Range: 3.125-100 ng/mL