

For Research Use Only

# LIPA Recombinant antibody, PBS Only (Capture)

Catalog Number: 85607-3-PBS



## Basic Information

Catalog Number:

85607-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3627

GenBank Accession Number:

BC012287

GeneID (NCBI):

3988

UNIPROT ID:

P38571

Full Name:

lipase A, lysosomal acid, cholesterol esterase

Calculated MW:

399 aa, 45 kDa

Purification Method:

Protein A purification

CloneNo.:

243025D11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01970-1, LIPA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85607-3-PBS. Detection antibody: 85607-2-PBS. Standard: Ag3627. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01970-2, LIPA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85607-3-PBS. Detection antibody: 85607-1-PBS. Standard: Ag3627. Range: 0.781-100 ng/mL.