

For Research Use Only

FLT3LG Recombinant antibody, PBS Only (Capture)

Catalog Number: 85656-5-PBS



Basic Information

Catalog Number:

85656-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001459.3

GeneID (NCBI):

2323

UNIPROT ID:

P49771-1

Full Name:

fms-related tyrosine kinase 3 ligand

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242884D2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02023-3, Human FLT3 Ligand Recombinant Matched Antibody Pair - PBS only. 85656-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2926. 85656-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



Serum of six individual healthy human donors was measured. The FLT3 Ligand concentration of detected samples was determined to be 179.5 pg/mL with a range of 89.9-235.2 pg/mL