

For Research Use Only

SF3B4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85663-2-PBS



Basic Information

Catalog Number: 85663-2-PBS	GenBank Accession Number: BC004273	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 10262	CloneNo.: 243098B6
Source: Rabbit	UNIPROT ID: Q15427	
Isotype: IgG	Full Name: splicing factor 3b, subunit 4, 49kDa	
Immunogen Catalog Number: AG0737	Calculated MW: 44 aa, 3 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02012-1, SF3B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85663-3-PBS. Detection antibody: 85663-2-PBS. Standard: Ag0737. Range: 0.781-100 ng/mL.