

For Research Use Only

SII-TFIIS,TCEA1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85664-3-PBS



Basic Information

Catalog Number: 85664-3-PBS	GenBank Accession Number: BC072460	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 6917	CloneNo.: 243155C9
Source: Rabbit	UNIPROT ID: P23193	
Isotype: IgG	Full Name: transcription elongation factor A (SII), 1	
Immunogen Catalog Number: AG12199	Calculated MW: 301 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2021-2, TCEA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85664-3-PBS. Detection antibody: 85664-2-PBS. Standard: Ag12199. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MPO2021-1, TCEA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85664-1-PBS. Detection antibody: 85664-3-PBS. Standard: Ag12199. Range: 0.781-50 ng/mL.