

For Research Use Only

APOA2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85667-1-PBS



Basic Information

Catalog Number:

85667-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001643.1

GeneID (NCBI):

336

UNIPROT ID:

P02652

Full Name:

apolipoprotein A-II

Calculated MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

243102C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02042-1, Apolipoprotein A II/APOA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85667-2-PBS. Detection antibody: 85667-1-PBS. Standard: Eg3085. Range: 0.781-100 ng/mL