

For Research Use Only

SMCR7,MID49 Recombinant antibody, PBS Only (Detector)

Catalog Number:85694-3-PBS



Basic Information

Catalog Number: 85694-3-PBS	GenBank Accession Number: BC014973	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 125170	CloneNo.: 243061C5
Source: Rabbit	UNIPROT ID: Q96C03	
Isotype: IgG	Full Name: Smith-Magenis syndrome chromosome region, candidate 7	
Immunogen Catalog Number: AG30119	Calculated MW: 454 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02025-1, SMCR7/MID49 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85694-4-PBS. Detection antibody: 85694-3-PBS. Standard: Ag30119. Range: 0.781-100 ng/mL