

For Research Use Only

A2BP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85732-4-PBS



Basic Information

Catalog Number:

85732-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18527

GenBank Accession Number:

BC113691

GeneID (NCBI):

54715

UNIPROT ID:

Q9NWB1

Full Name:

ataxin 2-binding protein 1

Calculated MW:

395 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

243149D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02075-2, A2BP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85732-3-PBS. Detection antibody: 85732-4-PBS. Standard: Ag18527. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP02075-1, A2BP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85732-2-PBS. Detection antibody: 85732-4-PBS. Standard: Ag18527. Range: 0.469-60 ng/mL