

For Research Use Only

C1S Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85736-5-PBS



Basic Information

Catalog Number:

85736-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5977

GenBank Accession Number:

BC056903

GeneID (NCBI):

716

UNIPROT ID:

P09871

Full Name:

complement component 1, s
subcomponent

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

250031G6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of eight individual healthy human donors was measured. The C1S concentration of detected samples was determined to be 610.10 ng/mL with a range of 102.92-1,188.99 ng/mL.



Sandwich ELISA standard curve of MP02077-3, Human C1S Recombinant Matched Antibody Pair - PBS only. 85736-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5977. 85736-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.