

For Research Use Only

# IL-12A/IL-12 p35 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85754-2-PBS



## Basic Information

Catalog Number:

85754-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC104982

GeneID (NCBI):

3592

UNIPROT ID:

P29459

Full Name:

interleukin 12A (natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35)

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

242887A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP02056-1, IL-12A/IL-12 p35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85754-3-PBS. Detection antibody: 85754-2-PBS. Standard: Eg3341. Range: 0.781-50 ng/mL.